



Sherer Zirconia

providing you

STRENGTH
&
Beauty



	FCZ Choice	FCZ Premium	FCZ Elite
Esthetics	☆☆☆	☆☆☆☆	☆☆☆☆☆
Strength	1100 - 1200 MPa 3Y	700 - 900 MPa 4Y	600 - 750 MPa 5Y
Recommended Usage	Posterior	Anterior or Posterior	Anterior or Posterior
Bridges	Long Span Bridges	Up to 3 - 4 Unit Bridges	Not recommended for bridges
Benefits	Superior strength allowing for maximum protection against bruxism	Increased Translucency for Optimal Esthetics	Multi-layer esthetics providing a natural transition from body to incisal

Call today to discuss the best zirconia options for your patients!



800-845-1116

www.shererdentallab.com